




Recombinant Mouse LIM and senescent cell antigen-like-containing domain protein 1 (Lims1)

Product Code	CSB-BP859129MO
Abbreviation	Recombinant Mouse Lims1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99JW4
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ANALASATCERCKGGFAPAEEKIVNSNGELYHEQCFVCAQCQFQQFPEGLFYEF EGRKYCEHDFQMLFAPCCHQCGEFIIGRVIKAMNNNSWHPECFRCDLCQEVLA DIGFVKNAGRHLCPCHNREKARGLGKYICQKCHAIIDEQPLIFKNDPYHPDHF NCANCGKELTADARELKGELYCLPCHDKMGVPICGACRRPIEGRVVNAMGKQ WHVEHFVCAKCEKPF LGHRHYERKGLAYCETHYNQLFGDVCFHCNRVIEGDV VSALNKAWCVSCFACSTCNTKLT LKNKFVEFDMKPVCKKCYEKFPLELKKRLK KLSETLGRK
Research Area	Cancer
Source	Baculovirus
Target Names	Lims1
Protein Names	Recommended name: LIM and senescent cell antigen-like-containing domain protein 1 Alternative name(s): Particularly interesting new Cys-His protein 1 Short name= PINCH-1
Expression Region	2-325aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	39.6 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.