




CDKL5 Antibody

Product Code	CSB-PA005085LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76039
Immunogen	Recombinant Human Cyclin-dependent kinase-like 5 protein (790-921AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase-like 5 (EC 2.7.11.22) (Serine/threonine-protein kinase 9), CDKL5, STK9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CDKL5

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, PC3 whole cell lysate, 293 whole cell lysate, K562 whole cell lysate, HepG2 whole cell lysate, Rat spleen tissue

All lanes: CDKL5 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 116, 108 kDa

Observed band size: 116 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.