



# FOXO4 Antibody, FITC conjugated

|                            |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA008839OC01HU                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                               |
| <b>Uniprot No.</b>         | P98177                                                                                      |
| <b>Immunogen</b>           | Peptide sequence from Human Forkhead box protein O4 protein (45-63AA)                       |
| <b>Raised In</b>           | Rabbit                                                                                      |
| <b>Species Reactivity</b>  | Human                                                                                       |
| <b>Form</b>                | Liquid                                                                                      |
| <b>Conjugate</b>           | FITC                                                                                        |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4            |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                   |
| <b>Isotype</b>             | IgG                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                  |
| <b>Alias</b>               | Forkhead box protein O4 (Fork head domain transcription factor AFX1), FOXO4, AFX AFX1 MLLT7 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                        |
| <b>Research Area</b>       | Cell Biology                                                                                |
| <b>Target Names</b>        | FOXO4                                                                                       |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                 |