



FOXO4 Antibody, Biotin conjugated

Product Code	CSB-PA008839OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98177
Immunogen	Peptide sequence from Human Forkhead box protein O4 protein (45-63AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein O4 (Fork head domain transcription factor AFX1), FOXO4, AFX AFX1 MLLT7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FOXO4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.