



FRYL Antibody, HRP conjugated

Product Code	CSB-PA009007LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94915
Immunogen	Recombinant Human Protein furry homolog-like protein (42-739AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein furry homolog-like (ALL1-fused gene from chromosome 4p12 protein), FRYL, AF4P12 KIAA0826
Immunogen Species	Homo sapiens (Human)
Target Names	FRYL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.