



SKOR2 Antibody, HRP conjugated

Product Code	CSB-PA009018PB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2VWA4
Immunogen	Peptide sequence from Human SKI family transcriptional corepressor 2 protein (19-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SKI family transcriptional corepressor 2 (Functional Smad-suppressing element on chromosome 18) (Fussel-18) (LBX1 corepressor 1-like protein) (Ladybird homeobox corepressor 1-like protein), SKOR2, CORL2 FUSSEL18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SKOR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.