



# SKOR2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA009018PD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q2VWA4
<b>Immunogen</b>	Peptide sequence from Human SKI family transcriptional corepressor 2 protein (19-36AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	SKI family transcriptional corepressor 2 (Functional Smad-suppressing element on chromosome 18) (Fussel-18) (LBX1 corepressor 1-like protein) (Ladybird homeobox corepressor 1-like protein), SKOR2, CORL2 FUSSEL18
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SKOR2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.