



GPLD1 Antibody, Biotin conjugated

Product Code	CSB-PA009721LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80108
Immunogen	Recombinant Human Phosphatidylinositol-glycan-specific phospholipase D protein (310-565AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol-glycan-specific phospholipase D (PI-G PLD) (EC 3.1.4.50) (Glycoprotein phospholipase D) (Glycosyl-phosphatidylinositol-specific phospholipase D) (GPI-PLD) (GPI-specific phospholipase D), GPLD1, PIGPLD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GPLD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.