



LPAR3 Antibody, HRP conjugated

Product Code	CSB-PA013049OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBY5
Immunogen	Peptide sequence from Human Lysophosphatidic acid receptor 3 protein (312-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 3 (LPA receptor 3) (LPA-3) (Lysophosphatidic acid receptor Edg-7), LPAR3, EDG7 LPA3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LPAR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.