



LPAR5 Antibody, HRP conjugated

Product Code	CSB-PA013051OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1C0
Immunogen	Peptide sequence from Human Lysophosphatidic acid receptor 5 protein (310-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 5 (LPA receptor 5) (LPA-5) (G-protein coupled receptor 92) (G-protein coupled receptor 93), LPAR5, GPR92 GPR93
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPAR5