



LPAR5 Antibody, Biotin conjugated

Product Code	CSB-PA013051OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1C0
Immunogen	Peptide sequence from Human Lysophosphatidic acid receptor 5 protein (310-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 5 (LPA receptor 5) (LPA-5) (G-protein coupled receptor 92) (G-protein coupled receptor 93), LPAR5, GPR92 GPR93
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPAR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.