



MST1R Antibody, Biotin conjugated

Product Code	CSB-PA015056LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04912
Immunogen	Recombinant Human Macrophage-stimulating protein receptor protein (1175-1400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage-stimulating protein receptor (MSP receptor) (EC 2.7.10.1) (CDw136) (Protein-tyrosine kinase 8) (p185-Ron) (CD antigen CD136) [Cleaved into: Macrophage-stimulating protein receptor alpha chain; Macrophage-stimulating protein receptor beta chain], MST1R, PTK8 RON
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MST1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.