



PRICKLE3 Antibody, HRP conjugated

Product Code	CSB-PA018678PB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43900
Immunogen	Peptide sequence from Human Prickle-like protein 3 protein (9-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prickle planar cell polarity protein 3 (LIM domain only protein 6) (LMO-6) (Prickle-like protein 3) (Pk3) (Triple LIM domain protein 6), PRICKLE3, LMO6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRICKLE3