



PRICKLE3 Antibody, Biotin conjugated

Product Code	CSB-PA018678PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43900
Immunogen	Peptide sequence from Human Prickle-like protein 3 protein (9-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prickle planar cell polarity protein 3 (LIM domain only protein 6) (LMO-6) (Prickle-like protein 3) (Pk3) (Triple LIM domain protein 6), PRICKLE3, LMO6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRICKLE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.