



PSMD14 Antibody, Biotin conjugated

Product Code	CSB-PA018904OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00487
Immunogen	Peptide sequence from Human 26S proteasome non-ATPase regulatory subunit 14 protein (264-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 14 (EC 3.4.19.-) (26S proteasome regulatory subunit RPN11) (26S proteasome-associated PAD1 homolog 1), PSMD14, POH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.