



PSMD4 Antibody, Biotin conjugated

Product Code	CSB-PA018908HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55036
Immunogen	Recombinant Human 26S proteasome non-ATPase regulatory subunit 4 protein (161-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 4 (26S proteasome regulatory subunit RPN10) (26S proteasome regulatory subunit S5A) (Antisecretory factor 1) (AF) (ASF) (Multiubiquitin chain-binding protein), PSMD4, MCB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.