



PTK7 Antibody, Biotin conjugated

Product Code	CSB-PA018996LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13308
Immunogen	Recombinant Human Inactive tyrosine-protein kinase 7 protein (772-1070AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive tyrosine-protein kinase 7 (Colon carcinoma kinase 4) (CCK-4) (Protein-tyrosine kinase 7) (Pseudo tyrosine kinase receptor 7) (Tyrosine-protein kinase-like 7), PTK7, CCK4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PTK7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.