



PTBP3 Antibody, FITC conjugated

Product Code	CSB-PA020060NC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95758
Immunogen	Peptide sequence from Human Polypyrimidine tract-binding protein 3 protein (37-55AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polypyrimidine tract-binding protein 3 (Regulator of differentiation 1) (Rod1), PTBP3, ROD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PTBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.