



PTBP3 Antibody, Biotin conjugated

Product Code	CSB-PA020060ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95758
Immunogen	Peptide sequence from Human Polypyrimidine tract-binding protein 3 protein (37-55AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polypyrimidine tract-binding protein 3 (Regulator of differentiation 1) (Rod1), PTBP3, ROD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PTBP3