



SNCA Antibody, Biotin conjugated

Product Code	CSB-PA021912OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37840
Immunogen	Peptide sequence from Human Alpha-synuclein protein (119-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-synuclein (Non-A beta component of AD amyloid) (Non-A4 component of amyloid precursor) (NACP), SNCA, NACP PARK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.