



TP73 Antibody, Biotin conjugated

Product Code	CSB-PA024092OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15350
Immunogen	Peptide sequence from Human Tumor protein p73 protein (537-555AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p73 (p53-like transcription factor) (p53-related protein), TP73, P73
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP73
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.