



UBE2M Antibody, HRP conjugated

Product Code	CSB-PA025468OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61081
Immunogen	Peptide sequence from Human NEDD8-conjugating enzyme Ubc12 protein (11-28AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NEDD8-conjugating enzyme Ubc12 (EC 2.3.2.-) (NEDD8 carrier protein) (Ubiquitin-conjugating enzyme E2 M), UBE2M, UBC12
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.