



# UBE2M Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA025468OC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P61081
<b>Immunogen</b>	Peptide sequence from Human NEDD8-conjugating enzyme Ubc12 protein (11-28AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	NEDD8-conjugating enzyme Ubc12 (EC 2.3.2.-) (NEDD8 carrier protein) (Ubiquitin-conjugating enzyme E2 M), UBE2M, UBC12
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	UBE2M
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.