



UBE2M Antibody, Biotin conjugated

Product Code	CSB-PA025468OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61081
Immunogen	Peptide sequence from Human NEDD8-conjugating enzyme Ubc12 protein (11-28AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NEDD8-conjugating enzyme Ubc12 (EC 2.3.2.-) (NEDD8 carrier protein) (Ubiquitin-conjugating enzyme E2 M), UBE2M, UBC12
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2M