



GLU1 Antibody, Biotin conjugated

Product Code	CSB-PA24649D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ZNZ7
Immunogen	Recombinant Arabidopsis thaliana Ferredoxin-dependent glutamate synthase 1, chloroplastic/mitochondrial protein (106-505AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ferredoxin-dependent glutamate synthase 1, chloroplastic/mitochondrial (EC 1.4.7.1) (Fd-GOGAT 1), GLU1, GLS1
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	GLU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.