



HPR Antibody, Biotin conjugated

Product Code	CSB-PA24659D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C9W5
Immunogen	Recombinant Arabidopsis thaliana Glycerate dehydrogenase HPR, peroxisomal protein (1-386AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerate dehydrogenase HPR, peroxisomal (GDH) (EC 1.1.1.29) (NADH-dependent hydroxypyruvate reductase 1) (AtHPR1) (HPR 1), HPR
Immunogen Species	Arabidopsis thaliana
Research Area	Cardiovascular
Target Names	HPR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.