




SUN1 Antibody

Product Code	CSB-PA24689A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FF75
Immunogen	Recombinant Arabidopsis thaliana SUN domain-containing protein 1 protein (2-105AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SUN domain-containing protein 1 (AtSUN1), SUN1
Immunogen Species	Arabidopsis thaliana
Research Area	Tags & Cell Markers
Target Names	SUN1

Image



Western Blot

Positive WB detected in: Spinacia oleracea leaf tissue

All lanes: SUN1 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.