



RB38 Antibody, Biotin conjugated

Product Code	CSB-PA24729D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6EMK7
Immunogen	Recombinant Chlamydomonas reinhardtii 38 kDa RNA-binding protein (1-382AA)
Raised In	Rabbit
Species Reactivity	Chlamydomonas reinhardtii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RB38
Immunogen Species	Chlamydomonas reinhardtii
Research Area	Others
Target Names	RB38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.