



TGFBR3L Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA27472B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | H3BV60 |
| Immunogen | Peptide sequence from Human Transforming growth factor-beta receptor type 3-like protein (242-258AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transforming growth factor-beta receptor type 3-like protein (TGF-beta receptor type-3-like protein) (TGFR-3L) (Transforming growth factor-beta receptor type III-like protein) (TGF-beta receptor type III-like protein), TGFBR3L |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | TGFBR3L |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |