



TGFBR3L Antibody, Biotin conjugated

Product Code	CSB-PA27472D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	H3BV60
Immunogen	Peptide sequence from Human Transforming growth factor-beta receptor type 3-like protein (242-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor-beta receptor type 3-like protein (TGF-beta receptor type-3-like protein) (TGFR-3L) (Transforming growth factor-beta receptor type III-like protein) (TGF-beta receptor type III-like protein), TGFBR3L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TGFBR3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.