



# PSG9 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA27499B0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q00887
<b>Immunogen</b>	Recombinant Human Pregnancy-specific beta-1-glycoprotein 9 protein (241-426AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pregnancy-specific beta-1-glycoprotein 9 (PS-beta-G-9) (PSBG-9) (Pregnancy-specific glycoprotein 9) (PS34) (Pregnancy-specific beta-1 glycoprotein B) (PS-beta-B) (Pregnancy-specific beta-1-glycoprotein 11) (PS-beta-G-11) (PSBG-11) (Pregnancy-specific glycoprotein 11) (Pregnancy-specific glycoprotein 7) (PSG7), PSG9, PSG11
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	PSG9
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.