




CHN48 Antibody

Product Code	CSB-PA362364LA01NHE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08252
Immunogen	Recombinant Nicotiana tabacum Endochitinase A protein (24-322AA)
Raised In	Rabbit
Species Reactivity	Nicotiana tabacum
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endochitinase A (CHN-A) (EC 3.2.1.14), CHN48
Immunogen Species	Nicotiana tabacum
Research Area	Others
Target Names	CHN48

Image



Western Blot

Positive WB detected in: Cucurbita moschata leaf tissue, Solanum tuberosum leaf tissue, Solanum lycopersicum leaf tissue

All lanes: CHN48 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36 kDa

Observed band size: 36 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.