



CHN48 Antibody, FITC conjugated

Product Code	CSB-PA362364LC01NHE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08252
Immunogen	Recombinant Nicotiana tabacum Endochitinase A protein (24-322AA)
Raised In	Rabbit
Species Reactivity	Nicotiana tabacum
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endochitinase A (CHN-A) (EC 3.2.1.14), CHN48
Immunogen Species	Nicotiana tabacum
Research Area	Others
Target Names	CHN48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.