



PETE Antibody, HRP conjugated

Product Code	CSB-PA365424HB01FKI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00289
Immunogen	Recombinant Spinacia oleracea Plastocyanin, chloroplastic protein (70-168AA)
Raised In	Rabbit
Species Reactivity	Spinacia oleracea
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plastocyanin, chloroplastic, PETE
Immunogen Species	Spinacia oleracea
Research Area	Others
Target Names	PETE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.