



rbcl Antibody, HRP conjugated

Product Code	CSB-PA517305LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O03042
Immunogen	Recombinant Arabidopsis thaliana Ribulose biphosphate carboxylase large chain protein (3-479AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribulose biphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	rbcl