



rbcl Antibody, FITC conjugated

Product Code	CSB-PA517305LC01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O03042
Immunogen	Recombinant Arabidopsis thaliana Ribulose biphosphate carboxylase large chain protein (3-479AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribulose bisphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	rbcl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.