



rbcl Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA517305LD01DOA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O03042 |
| Immunogen | Recombinant Arabidopsis thaliana Ribulose biphosphate carboxylase large chain protein (3-479AA) |
| Raised In | Rabbit |
| Species Reactivity | Arabidopsis thaliana |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ribulose biphosphate carboxylase large chain (RuBisCO large subunit) (EC 4.1.1.39), rbcl |
| Immunogen Species | Arabidopsis thaliana |
| Research Area | Others |
| Target Names | rbcl |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |