



# PSMD6 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA614875PB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15008
<b>Immunogen</b>	Peptide sequence from Human 26S proteasome non-ATPase regulatory subunit 6 protein (355-372AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	26S proteasome non-ATPase regulatory subunit 6 (26S proteasome regulatory subunit RPN7) (26S proteasome regulatory subunit S10) (Breast cancer-associated protein SGA-113M) (Phosphonoformate immuno-associated protein 4) (Proteasome regulatory particle subunit p44S10) (p42A), PSMD6, KIAA0107 PFAAP4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	PSMD6
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.