



PSMD6 Antibody, FITC conjugated

Product Code	CSB-PA614875PC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15008
Immunogen	Peptide sequence from Human 26S proteasome non-ATPase regulatory subunit 6 protein (355-372AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 6 (26S proteasome regulatory subunit RPN7) (26S proteasome regulatory subunit S10) (Breast cancer-associated protein SGA-113M) (Phosphonoformate immuno-associated protein 4) (Proteasome regulatory particle subunit p44S10) (p42A), PSMD6, KIAA0107 PFAAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.