



AOX1A Antibody, Biotin conjugated

Product Code	CSB-PA655106LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39219
Immunogen	Recombinant Arabidopsis thaliana Ubiquinol oxidase 1a, mitochondrial protein (63-178AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquinol oxidase 1a, mitochondrial (EC 1.10.3.11) (Alternative oxidase 1a), AOX1A, AOX1 HSR3
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	AOX1A