



PPM1J Antibody, HRP conjugated

Product Code	CSB-PA686382LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JR12
Immunogen	Recombinant Human Protein phosphatase 1J protein (294-505AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1J (EC 3.1.3.16) (Protein phosphatase 2C isoform zeta) (PP2C-zeta), PPM1J, PPP2CZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPM1J
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.