



CESA4 Antibody, HRP conjugated

Product Code	CSB-PA771138PB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q84JA6
Immunogen	Peptide sequence from Arabidopsis thaliana Cellulose synthase A catalytic subunit 4 [UDP-forming] protein (924-937AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellulose synthase A catalytic subunit 4 [UDP-forming] (AtCesA4) (EC 2.4.1.12) (Protein IRREGULAR XYLEM 5) (AtIRX5), CESA4, IRX5
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	CESA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.