



RTN4RL2 Antibody, Biotin conjugated

Product Code	CSB-PA774805ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UN3
Immunogen	Peptide sequence from Human Reticulon-4 receptor-like 2 protein (355-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4 receptor-like 2 (Nogo receptor-like 3) (Nogo-66 receptor homolog 1) (Nogo-66 receptor-related protein 2) (NgR2), RTN4RL2, NGRH1 NGRL3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RTN4RL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.