



RTN4RL2 Antibody, Biotin conjugated

Product Code	CSB-PA774805ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UN3
Immunogen	Peptide sequence from Human Reticulon-4 receptor-like 2 protein (355-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reticulon-4 receptor-like 2 (Nogo receptor-like 3) (Nogo-66 receptor homolog 1) (Nogo-66 receptor-related protein 2) (NgR2), RTN4RL2, NGRH1 NGRL3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RTN4RL2