




AMOTL2 Antibody

Product Code	CSB-PA809099NA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2J4
Immunogen	Peptide sequence from Human Angiotenin-like protein 2 protein (718-737AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiotenin-like protein 2 (Leman coiled-coil protein) (LCCP), AMOTL2, KIAA0989
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AMOTL2

Image



Western Blot

Positive WB detected in: 293T whole cell lysate, U251 whole cell lysate, HepG2 whole cell lysate, 293 whole cell lysate, A549 whole cell lysate, MCF-7 whole cell lysate, PC3 whole cell lysate
All lanes: AMOTL2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 86, 92 kDa
Observed band size: 86 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.