




DEG15 Antibody

Product Code	CSB-PA819818LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VZD4
Immunogen	Recombinant Arabidopsis thaliana Glyoxysomal processing protease, glyoxysomal protein (572-709AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyoxysomal processing protease, glyoxysomal (AtDEG15) (EC 3.4.21.-) (DEG-protease), DEG15, GPP
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	DEG15

Image



Western Blot

Positive WB detected in: Spinacia oleracea leaf tissue

All lanes: DEG15 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 77 kDa

Observed band size: 77 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.