



CAT1 Antibody, HRP conjugated

Product Code	CSB-PA822149LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96528
Immunogen	Recombinant Arabidopsis thaliana Catalase-1 protein (1-492AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Catalase-1 (EC 1.11.1.6), CAT1
Immunogen Species	Arabidopsis thaliana
Research Area	Neuroscience
Target Names	CAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.