



DNAJA3 Antibody, Biotin conjugated

Product Code	CSB-PA836207LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EY1
Immunogen	Recombinant Human DnaJ homolog subfamily A member 3, mitochondrial protein (86-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily A member 3, mitochondrial (DnaJ protein Tid-1) (hTid-1) (Hepatocellular carcinoma-associated antigen 57) (Tumorous imaginal discs protein Tid56 homolog), DNAJA3, HCA57 TID1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJA3