



GLO1 Antibody, HRP conjugated

Product Code	CSB-PA862860LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9LRR9
Immunogen	Recombinant Arabidopsis thaliana (S)-2-hydroxy-acid oxidase GLO1 protein (1-367AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	(S)-2-hydroxy-acid oxidase GLO1 (EC 1.1.3.15) (Glycolate oxidase 1) (AtGLO1) (GOX 1) (Short chain alpha-hydroxy acid oxidase GLO1), GLO1
Immunogen Species	Arabidopsis thaliana
Research Area	Signal Transduction
Target Names	GLO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.