



# GLO1 Antibody, FITC conjugated

|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA862860LC01DOA                                                                                                                      |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                            |
| <b>Uniprot No.</b>         | Q9LRR9                                                                                                                                   |
| <b>Immunogen</b>           | Recombinant Arabidopsis thaliana (S)-2-hydroxy-acid oxidase GLO1 protein (1-367AA)                                                       |
| <b>Raised In</b>           | Rabbit                                                                                                                                   |
| <b>Species Reactivity</b>  | Arabidopsis thaliana                                                                                                                     |
| <b>Form</b>                | Liquid                                                                                                                                   |
| <b>Conjugate</b>           | FITC                                                                                                                                     |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                         |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                               |
| <b>Alias</b>               | (S)-2-hydroxy-acid oxidase GLO1 (EC 1.1.3.15) (Glycolate oxidase 1) (AtGLO1) (GOX 1) (Short chain alpha-hydroxy acid oxidase GLO1), GLO1 |
| <b>Immunogen Species</b>   | Arabidopsis thaliana                                                                                                                     |
| <b>Research Area</b>       | Signal Transduction                                                                                                                      |
| <b>Target Names</b>        | GLO1                                                                                                                                     |